




Sequencing technologies



Illumina: Sequencing Platforms

<https://www.illumina.com/systems/sequencing-platforms.html>

Benchtop



iSeq 100 System



MiniSeq System



MiSeq Series



NextSeq Series

Run Time	9–17.5 hours	4–24 hours	4–55 hours	12–30 hours
Maximum Output	1.2 Gb	7.5 Gb	15 Gb	120 Gb
Maximum Reads Per Run	4 million	25 million	25 million [†]	400 million
Maximum Read Length	2 × 150 bp	2 × 150 bp	2 × 300 bp	2 × 150 bp



NextSeq Series



HiSeq Series



HiSeq X Series[‡]



NovaSeq 6000 System


Production-Scale

Run Time	12–30 hours	< 1–3.5 days (HiSeq 3000/HiSeq 4000) 7 hours–6 days (HiSeq 2500)	< 3 days	16–36 hours (Dual S2 flow cells) 44 hours (Dual S2 flow cells)
Maximum Output	120 Gb	1500 Gb	1800 Gb	6000 Gb
Maximum Reads Per Run	400 million	5 billion	6 billion	20 billion
Maximum Read Length	2 × 150 bp	2 × 150 bp	2 × 150 bp	2 × 150 bp

Flow cell

- ▶ A glass surface coated with an arrangement of paired oligos that are complementary to the adapters added to your template molecules.
- ▶ A flow cell can have varying numbers of lanes, depending on the sequencing machine.

http://training.bioinformatics.ucdavis.edu/docs/2014/09/september-2014-workshop/Monday_JF_HTS_lecture.html



HiSeq 2500 (2 lane)




HiSeq 3000/4000



NextSeq 500

TTTATGATACGGCGACCGAGAUCTACAC-3'

TTTCAAGCAGAAGACGCATACGAGGxxAT-3'



Illumina: flow cell

Introduction to Sequencing by Synthesis

Video courtesy of [Illumina YouTube channel](#)

Introduction to Sequencing by Synthesis


Video courtesy of [Illumina YouTube channel](#)

Number of clusters \sim Number of reads


Number of sequencing cycles \sim Length of reads

Single end or Paired end?

Depends on your goals, paired-end reads are better for reads that map to multiple locations, for assemblies and for splice isoform differentiation.




Options for sequencing




- ✓ SE - Single end dataset => Only Read1

Options for sequencing



- ✓ SE - Single end dataset => Only Read1
- ✓ PE - Paired-end dataset => Read1 + Read2
 - can be 2 separate FastQ files or just one with interleaved pairs

Options for sequencing



- ✓ SE - Single end dataset => Only Read1
- ✓ PE - Paired-end dataset => Read1 + Read2
 - can be 2 separate FastQ files or just one with interleaved pair
- ✓ Fragment length: ~300-500bp
- ✓ Read length: 50bp - 250bp, depends on the sequencer (HiSeq2500, MiSeq, NextSeq)

Options for sequencing

Single end or Paired end?

For differential gene expression, which one you pick depends on-

- If you are specifically interested in **isoform-level differences**
- The abundance of **paralogous genes** in your system of interest
- Your **budget**, paired-end data is usually 2x more expensive


Multiplexing (with barcodes and indices)

Multiplexing (with barcodes and indices)

- ✦ Charges for sequencing are usually per lane of the flow cell, and usually you don't need one lane per sample

Multiplexing (with barcodes and indices)

- ✦ Charges for sequencing are usually per lane of the flow cell, and usually you don't need one lane per sample
- ✦ Multiplexing allows you to sequence multiple samples per lane with addition of indices (within the Illumina adapter) or special barcodes (outside the Illumina adapter).



Pacific Biosciences: <http://www.pacb.com/>

Oxford Nanopore (MinION): <https://nanoporetech.com/>

10X Genomics: <https://www.10xgenomics.com/>

Other Sequencing Platforms

	Advantages	Disadvantages
<u>Pacific Biosciences</u>	Iso-Seq protocol for transcripts up to 10Kb, high base calling accuracy	High cost, large machines
<u>Oxford Nanopore</u>	Accurate quantitative data for short transcripts (< 700bp), portable, high yield	High errors rate affects assembling de novo transcripts, higher amount of cDNA input
<u>10X Genomics</u>	Low cost (integrated with short-read technology), barcoding for accurate isoform detection, low error rates	Extra preparation step (barcode), extra computational step

Transcriptomics with long read technologies

	Advantages	Disadvantages
<u>Pacific Biosciences</u>	Iso-Seq protocol for transcripts up to 10Kb, high base calling accuracy	High cost, large machines
<u>Oxford Nanopore</u>	Direct RNA-seq kit, portable, high yield, PCR-free	High errors rate affects assembling de novo transcripts, higher amount of cDNA input
<u>10X Genomics</u>	Low cost (integrated with short-read technology), barcoding for accurate isoform detection, low error rates	Extra preparation step (barcode), extra computational step

Transcriptomics with long read technologies

These materials have been developed by members of the teaching team at the Harvard Chan Bioinformatics Core (HBC). These are open access materials distributed under the terms of the Creative Commons Attribution license (CC BY 4.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

